

# Hyperspectral imaging facilitates early detection of intraoperative ischaemia during neurovascular surgery

Oscar J MacCormac, MD<sup>1,2</sup>; Conor C Horgan, PhD<sup>3</sup>; Zhonghao Wang, MSc<sup>1</sup>; Matthew Elliot, MD<sup>1,2</sup>; Dale J Waterhouse, PhD<sup>3</sup>; Sebastien Ourselin, PhD<sup>1,3</sup>; Michael Ebner, PhD<sup>3</sup>, Christos Tolias, MD<sup>2</sup>; Ahilan Kailaya-Vasan, MD<sup>2</sup>; Tom Vercauteren, PhD<sup>2,3</sup>; Jonathan Shapey, MD, PhD<sup>1,2,3</sup> <sup>1</sup>School of Biomedical Engineering and Imaging Sciences, King's College London; <sup>2</sup>King's College Hospitals Dept. of Neurosurgery; <sup>3</sup>Hypervision Surgical Ltd.

B

### Introduction

Hyperspectral imaging (HSI) is an advanced, wide field optical imaging technique that enables the visualization of light reflectance beyond that of the naked eye.

In tissue, the refined spectral information from HSI permits real-time extraction of biological and physiological parameters, one of which is tissue oxygenation (StO<sub>2</sub>).

Intra-operative stroke occurs in up to 6.4% of cranial neurosurgical procedures, and thus early detection of intra-operative ischaemia may prove key to improving patient outcomes.








Figure 1. Snapshot mosaic hyperspectral camera mounted to Zeiss Kinevo-900 neurosurgical operating microscope

#### **Methods and Materials**

A snapshot mosaic HSI camera was coupled to a Zeiss Kinevo-900 (Fig. 1) and connected to a dedicated camera control unit.

HSI video recordings were acquired during eight aneurysm surgeries and eight arteriovenous malformation (AVM) extirpations.

Recordings manually annotated by a neurosurgeon (Fig 2.), facilitating extraction of StO<sub>2</sub> information from each annotated region of interest (RoI).

#### Results

Using our system, real-time HSI recordings were successfully and safely acquired during all 16 cases.

StO<sub>2</sub> information was extracted from tissue structures across cases.

HSI detected a decrease in StO<sub>2</sub> (Fig 3.) over a region of cortex 8 minutes before electrophysiological change was demonstrated during one MCA aneurysm clipping.






Figure 2. (A) HSI white light image during neurovascular surgery (B) Neurosurgeon annotated HSI recording from which StO<sub>2</sub> information derived




Figure 3. HSI derived StO<sub>2</sub> values for capillary and cortex at the beginning of surgery demonstrating healthy values (A) and ischaemic values (B) during inadvertent vessel occlusion

## **Conclusions**

Our HSI system has demonstrated capability in providing enhanced tissue information in real-time during surgery.

HSI is capable of detecting ischaemic changes much earlier than current methods, paving the way for earlier intra-operative decision making and safer surgery in the future.

